

Gel LECTRO for ultrasound - GUSLT5G-15-5RB



Description

Coupling gel for ultrasound transducers. High viscosity product with technical / chemical suitable for examinations:

- Cardiology, internist, obstetrics, gynecology
- Ultrasound, physiokinesitherapy, echo encephalography, Doppler.

Technical features

Code:	GUSLT5G-15-5RB
Name and description:	Gel LECTRO for ultrasound
Packaging:	5-liter rigid bag, Box of 4 bags
Composition:	Water, propylene glycol, carbomer, diazolidinyl urea, methylparaben, propylparaben, disodium EDTA
Color:	Transparent blue
Feature:	Manufactured with raw materials of pharmaceutical purity. pH=7. Water-soluble, non-greasy and non-staining, non-irritating, odourless, does not damage or corrode the probes, non-flammable.
Storing:	No problem of pollution or temperature. The stability of the product is greater than 30 months. Do not expose to strong UV radiation and do not place in direct contact with sources of heat.
Toxicity:	There are no cases of intolerance to the product when used according to the modes of operation described below. In case of accidental oral ingestion, the toxicity limit is 5 g / kg of body weight. It contains no perfumes, essences, formaldehyde.
Manufactured by:	TELIC S.A. - Spain

Mode of use

Screw to the mouth of the tank the dispenser supplied (1 dispenser each 4 tanks). Transfer the contents of the tank to fill the 250ml bottle supplied. Thoroughly and vigorously clean the patient's skin. Using the bottle, apply a small amount of gel on the affected area, and then use the probe. After examination eliminate any residue with paper towel. Thoroughly clean the probe after use.

Complying with norms

Medical device class I, Directive 93/42/EEC and subsequent amendments.
Italian Medical Devices Datastandard Registration: 48672/R; CND: V 9099