

Codice / Reference	Si raccorda con / It fits with	INT Ø 1	EXT Ø 1	INT Ø 2	EXT Ø 2
TSVBM001	Sensormedics - Cortex	27,0	30,0	27,0	34,5
TSVBM002	Micromedical - Vitalograph - MIR	27,0	30,0	30,5	30,0
TSVBM004	Cosmed - Customed	27,0	30,0	29,5	32,8
TSVBM005	Koko	27,0	30,0	45,5	48,0
TSVBM006	Medical Graphics	27,2	30,0	34,3	36,3

Disponibili con e senza boccaglio / Available with and without mouthpiece



Proteggete i vostri pazienti e gli strumenti con i filtri VBMAX. La bassa resistenza permette risultati molto accurati anche ai bassi flussi. I filtri VBMAX trattengono il 99,99% di batteri e virus e rispettano le raccomandazioni ATS circa la resistenza ai flussi nella strumentazione per test polmonari.

*Protect your patient and instruments with our VBMAX Filters! VBMAX provides the highest protection for the patient as well as for your equipment. Its low resistance allows an accurate testing. VBMAX filters are very effective and filters more than 99,99% of bacteria and viruses. At the same time, VBMAX easily surpasses the ATS recommendations about flow resistance.*

Confezionati singolarmente / Singolarly wrapped

## accessori per spirometria / spirometry accessories

### SIRINGA DI CALIBRAZIONE CALIBRATION SIRYNGE

Siringa di calibrazione per tutti gli spirometri sul mercato, ottimo rapporto qualità prezzo. Volume: 3 l

Calibration syringe for all the spirometers in the market, great value for the money. Volume: 3 l



*boccagli e filtri per spirometria page 7/8*



### DISTANZIATORE PER INALAZIONI LEVERHALER LEVERHALER METERED-DOSE INHALER

Tra i distanziatori presenti in commercio:

- Il più facile da usare
- Il più semplice da pulire
- Il più efficace
- Ha il prezzo maggiormente competitivo

Among the inhaler on the market:

- The easiest to use
- The easiest to clean
- The most effective
- The most competitive price

