

Codice /Reference	Si raccorda con / It fits with	INT Ø	EXT Ø
TS5218	Cosmed	28,0	30,5
TS3962	Micromedical - Vitalograph - MIR	28,0	30,0
TS3963	Fukuda Sangyo - Medical Graphics - Biomedin	25,3	27,0
TS52422	Customed	29,2	32,0
TS8182	Sensormedics - Microcard	30,3	33,0
TS8248	Siebel - BTL	26,2	28,0
TS8551	Pediatric standard - Smokerlyzer	18,5	22,5
TS9190	Fukuda Denshi - Welch Allyn	24,2	26,0
TS8822	FIDIM V/B Filters	25,4	26,8



DISPONIBILI CON E SENZA BOCCAGLIO / Available with and without mouthpiece
I PRIMI BOCCAGLI CONFEZIONATI SINGOLARMENTE / The first single wrapped mouthpieces

confezionamento:



dispenser 100 pz
dispenser 100 pcs



scatola 500 pz
box 500 pcs



pallet 45 scatole
pallet 45 boxes



siringa calibrazione spirometri ρ 8
calibration syringe ρ 8



richiedi la valigetta delle
campionature spirometria
ask for our spirometry
sample box

accessori per spirometria / accessories for spirometry

PINZETTE STRINGINASO **NOSECLIPS**

Pinzette stringinaso per esami spirometrici, confortevoli e delicate. Realizzate con morbidi materiali biocompatibili, assicurano la perfetta chiusura delle vie nasali durante ogni tipo di test respiratorio.



CONFEZIONAMENTO / PACKAGING
scatola/box da 10 pz/pcs

TURBINE MONOUSO FLOWMIR **FLOWMIR DISPOSABLE TURBINE**

FlowMir si adatta ai nuovi modelli di spirometri MIR, è già calibrata ed è testata in fabbrica con attrezzature sofisticate.



FlowMir disposable turbine fits all the MIR spirometers, calibrated and single packed.

CONFEZIONAMENTO / PACKAGING
100 pz/pcs