



"EUROCATCH" ECG Clamp Electrodes



DESCRIPTION

"EUROCATCH" LUMED Reusable Peripheral ECG Clamp Electrodes have features that make them reliable during signal recording with all ECG devices available on the market. The electrode metallic part is made of silver coated material in order to allow a safe and clear signal transmission. LUMED clamp electrodes, thanks to screw and clip connections, can be used with all ECG devices and monitors with patient cable that has banana pin (4mm) or snap. The product is latex and nickel free and the plastic clamp material is hypoallergenic. They are available in set of 4 clamp electrodes, as showed in the picture, and in the following versions: adult (EL0004) and paediatric (EL0005).

TECHNICAL FEATURES

REFERENCE	EL0005 <i>paediatric</i>	EL0004 <i>adult</i>
DESCRIPTION	Reusable peripheral ECG clamp electrodes	
ELECTRODE MATERIAL-COATING	METALLIC coated with SILVER/SILVER CHLORIDE (SSC)	
CONNECTION	UNIVERSAL - hole with screw (4 mm) and clip	
PLASTIC CLAMP MATERIAL	ABS NOVODUR BAYER	DURETHAN B30S (NYLON PA6)
COLOUR	(Yellow, green, red and black)	
CLAMP DIMENSIONS	20 x 85 mm	30 x 130 mm
COMPLIANCE TO STANDARDS	Directive 93/42/EEC and amendments. Class I medical device *Non-Active Medical Device*	
MANUFACTURER	LUMED srl Via Staffora, 18/9 - 20090 Opera (MI) - ITALY	
SALE PACKAGING	4 pcs in polyethylene pouch	
STORAGE	Temperature: 0°C + +50°C - Relative humidity: 20% + 80%	
PRODUCT SHELF LIFE	5 years from production date printed on packaging	