

Crocodile / Snap Adapters



DESCRIPTION

Thanks to LUMED crocodile adapter (ET101), you will be able to keep your patient cable with 4mm banana plug and use it with disposable electrodes that need a crocodile clip to be connected.

Thanks to LUMED snap adapter (ET102), you will be able to keep your patient cable with 4mm banana plug and use it with disposable electrodes that need a snap plug to be connected.

TECHINICAL FEATURES

REFERENCE	ET101	ET102
DESCRIPTION	Reusable crocodile / snap adapter for disposable electrodes	
CONNECTION TO CABLE	∅ 4mm socket	
CONNECTION SYSTEM	Crocodile	Female Snap
COLOUR	Red	White
MATERIALS	Plastic: PVC, PE - Metallic: copper, brass	
COMPLIANCE TO STANDARDS	Directive 93/42/EEC and amendments. Class I medical device <i>*Non-Active Medical Device*</i>	
MANUFACTURER	LUMED srl <i>Via Staffora, 18/9 - 20090 Opera (MI) - ITALY</i>	
SALE PACKAGING	10 same color pcs in zip-lock bag	
STORAGING	Temperature: 0°C ÷ +50°C – Relative humidity 20% ÷ 80%	
PRODUCT SHELF LIFE	5 years from production date if correctly store	