



# EUROBIKE

IS MADE IN LUMED®



 **LUMED®**  
[www.lumed.com](http://www.lumed.com)



**EUROBIKE** È UN PRODOTTO LUMED® / IS MADE IN LUMED®

Gli ergometri di precisione **EUROBIKE 320** ed **EUROBIKE 420** sono prodotti all'avanguardia e di elevate prestazioni. Il dispositivo risponde ai più elevati standard di qualità nei test per la esatta valutazione dello sforzo fisico nel campo della diagnostica della funzionalità cardiovascolare e polmonare.

Particolarmente indicato per la riabilitazione, grazie anche alla notevole accessibilità.

## CARATTERISTICHE

- Accesso facilitato per ogni tipo di paziente
- Robusta struttura in acciaio,
- Unità meccanica compatta
- Regolazione agevole del manubrio e della sella
- Intervallo di carico da 1 a 999 W
- Assoluta silenziosità grazie al sistema Poly-V®
- Misuratore di pressione arteriosa integrato (solo nel modello **EUROBIKE 420**)

## OPZIONI DISPONIBILI

- Configurazione pediatrica
- Frequenzimetro Polar
- Saturimetro integrato

Tutti i marchi registrati o depositati sono proprietà dei rispettivi titolari. I prodotti potrebbero differire da quanto illustrato. Il produttore si riserva il diritto di modificare i prodotti qui descritti senza darne preavviso.

Precision ergometers **EUROBIKE 320** and **420** are high performance enhanced devices, meeting the highest quality standards for cardio-pulmonary examinations. Especially suitable for Stress-Test and CPET, its easy accessibility makes **EUROBIKE** the best choice for Cardiac- and Physical Rehabilitation.

## MAIN FEATURES

- Easy access for any kind of patient
- Strong steel structure
- Compact mechanical unit
- Wide saddle and handlebar setting range
- Workload ranging from 1 to 999 W
- Noise-free thanx to the Poly-V® system
- Built-in automatic NIBP monitor (only for **EUROBIKE 420**)

## AVAILABLE OPTIONS

- Pediatric configuration
- Hear-Rate monitor
- SpO2 monitor

Every Trademark or Registered Mark is property of the respective owner. The actual products may differ from here shown. The Manufacturer has the right to modify the products without prior notice.



SEDE DI OPERA [MI]  
VIA STAFFORA 18/9 • 20090, OPERA [MI] • ITALY  
TEL: 02 57606751 • FAX: 02 57606752

SEDE DI FORLÌ [FC]  
VIA SENIO 36/40 • 47121, FORLÌ [FC] • ITALY  
TEL: 0543 702380 • FAX: 0543 707518