

Crocodile / Snap Adapters



DESCRIPTION

Thanks to LUMED crocodile adapter (ET101), you will be able to keep your patient cable with 4mm banana plug and use it with disposable electrodes that need a crocodile clip to be connected.

Thanks to LUMED snap adapter (ET102), you will be able to keep your patient cable with 4mm banana plug and use it with disposable electrodes that need a snap plug to be connected.

TECHINICAL FEATURES

REFERENCE	ET101	ET102
DESCRIPTION	Reusable crocodile / snap adapter for disposable electrodes	
COMPLIANCE TO STANDARDS	Directive 93/42/EEC and amendments. Class I medical device *Non-Active Medical Device*	
MANUFACTURER	LUMED srl Via Staffora, 18/9 - 20090 Opera (MI) - ITALY	