



"EUROSPIRO" Disposable Nose Clips



DESCRIPTION

"EUROSPIRO" Disposable Nose clips are intended for spirometry tests in order to reduce the risk of introducing a technical error. LUMED rubber- and foam-padded nose clips help prevent patients from inhaling or exhaling through the nose during spirometry.

LUMED nose clips are made of white thermoplastic, latex-free and CE marked. The tension in the nose clip body keeps it firmly in place while the rubber and foam pads are designed for the best patient comfort.

The nose clips are disposable and intended for single-patient use only.

TECHNICAL FEATURES

REFERENCE	TSTRA001/100	TSTRA001/10	TSTRA006/100	TSTRA006/10
DESCRIPTION	Disposable Rubber- and Foam-padded Nose Clips			
NOSE CLIP BODY MATERIAL	Polyacetate	Polyacetate	Polyacetate	Polyacetate
NOSE CLIP PADS MATERIAL	Polyvinyl	Polyvinyl	Ethylene-vinyl acetate	Ethylene-vinyl acetate
COLOUR	White			
DIMENSIONS	70 x 40 mm			
COMPLIANCE TO STANDARDS	<p>Directive 93/42/EEC and amendments. Class I medical device Italian Medical Device Data standard Registration: 187604/R; CND: R9099; GMDN: 10907; CE *Non-Active Medical Device*</p>			
MANUFACTURER	<p>LUMED srl Via Staffora, 18/9 - 20090 Opera (MI) - ITALY</p>			
SALE PACKAGING	100 pcs in polyethylene	10 pcs in cardboard box	100 pcs in polyethylene	10 pcs in cardboard box
STORAGE	Temperature: 0°C +50°C - Relative humidity: 20% + 80%			