


EE100CUBI-200 PATIENT CABLE FOR RESTING ECG EE100CUSI-200 PATIENT CABLE FOR STRESS-TESTING



NB the product might differ from the pictures here shown. The manufacturer can modify the product without prior notice.

Description

One-piece patient cable for Resting ECG and Stress-Testing (12-leads). It is completely compatible with the original ones.

The plastic connector towards the device is conformal to the original. The terminations have 4mm Banana plugs or Snaps.

The Compatibility List could be extended to further devices, according to the specific features.

Every trademark is propriety of the respective owner.

Technical Features

Length:	✓ Trunk: 200 cm ✓ Leadwires: chest 75 cm, plimbseriferiche 100 cm
Terminations:	✓ EE100CUBI-200: 10 terminali a banana 4 mm ✓ EE100CUSI-200: 10 terminali a snap
Yoke:	pressofuso, compatto
Connector:	tipo Siemens standard 15 pin, in plastica e con clip di fissaggio
Colour code:	IEC
Resistors:	4.7 kΩ con protezione contro le scariche di defibrillazione
Protection against defibrillation:	ASPEL AsCard 3 (1992 – 1993). BOSCH 603D, 603ER, 603S, 606. CAMBRIDGE CM6000, CM6400. CARDIOLINE Delta 300, 360. CARDIO CONTROL Cardio Perfect. CARDIOLEX EC Sense. CUSTO.MED card M, cardio 130. DEGO Cardiomed 3. HELIGE EK36, EK41, EK43, EK56, EK403, EK413, EK512, EK512 E, EK512 P. MARQUETTE HELIGE CardioSmart, CardioSys, Microsmart, MAC 500, MAC 1000, MAC 1200, Corina. HORMANN VON BERG Bioset 3500, 3600, 3700, 6000, 8000. SCHWARZER Cardos 6. SIEMENS Cardiostat 1, 3T, 703, 710, 31, 31S; Cardirex 3T, 4, 6T, 62; Megacard Elema A/B, R, R/E; Cathcor, Cardina C, Miracard 1; Sirecard 440, 440S, 460; Mingodoc, Mingocard 3, 4, 7; Mingophon 3, 7, 310, 410, 420, 710, 720; Mingograph 4, 62, 81, 82, 740; Mingolog 7; WELCH ALLYN Cardio Perfect. ZIMMER Cardio 1, Cardio 3 ...

