


## EE100FRI-200 PATIENT CABLE FOR RESTING ECG EE100FRSI-200 PATIENT CABLE FOR STRESS-TESTING



NB the product is subject to changes without prior notice.

### Description

One-piece patient cable for Resting ECG and Stress-Testing ECG (12-leads). It is completely compatible with the original ones.

The plastic connector towards the device is conformal to the original. The terminations have 4mm Banana or Snap plugs.

Made for devices from Fukuda Denshi.

The Compatibility List could be extended to further devices, according to the specific features.

Every trademark is propriety of the respective owner.

### Technical Features

|                       |   |
|-----------------------|---|
| <b>Length:</b>        | <ul style="list-style-type: none"> <li>✓ Trunk: 200 cm</li> <li>✓ Leadwires: chest 75 cm, limbs 100 cm</li> </ul>   |
| <b>Terminations:</b>  | <ul style="list-style-type: none"> <li>✓ EE100FRI-200: 4 mm Banana</li> <li>✓ EE100FRSI-200: Snap</li> </ul>  |
| <b>Yoke:</b>          | compact, sealed   |
| <b>Connector:</b>     | Fukuda style with clamps, made in plastic.  |
| <b>Colour code:</b>   | IEC   |
| <b>Resistors:</b>     | 10 kΩ with defi protection  |
| <b>Compatibility:</b> | <p>FUKUDA DENSHI Autocardiner: FCP11, FCP2101A, FCP220A, FCP4101A, FX101, MCP302.</p> <p>FUKUDA DENSHI Cardimax: FX2111, FCP2155, FCP7101, FX3010, FX4010, FX7102, FX7202, FX7302, FX7402, FX7542, FX8222, FX8322, FX8328.</p> <p>SIEMENS Cardiostat 1, 11, 31, 31S ...</p> |

