


EE100NRI-200 PATIENT CABLE FOR REST-ECG EE100NRSI-200 PATIENT CABLE FOR STRESS-TEST



NB the product is subject to changes without prior notice.

Description

One-piece patient cable for Rest-ECG (12-leads diagnostic ECG) and Stress-Test. It is completely compatible with the original ones.

The plastic connector towards the device is conformal to the original. The terminations have 4mm Banana plugs.

Made for devices from the following Brands: Nihon-Kohden (connector with screws. The Compatibility List could be extended to further devices, according to the specific features.

Every trademark is propriety of the respective owner.

Technical Features

Length:	<ul style="list-style-type: none"> ✓ Main cable: 200 cm ✓ Leadwires: chest 70 cm, limbs 110 cm
Terminations:	<ul style="list-style-type: none"> ✓ EE100NRI-200: 4 mm Banana ✓ EE100NRSI-200: Snap
Yoke:	compact, sealed
Connector:	D-sub 15 pin, made from plastic
Colour code:	IEC
Resistors:	10 k Ω
Protection against defibrillation:	present
Compatibility:	MEDSET Mobile Resting-ECG (Nihon-Kohden 1250, 1350, 1500, 1550). NIHON KHODEN ECG-1150 cardiofax C, 1250 cardiofax S, 1350 cardiofax M, 1500/-1550 cardiofax V, 2350 cardiofax M, PEA-1110. ECG-6511, 8110, 8240, 8270, 8330, 8340, 8350, 8370, 84xx, 882x, 883x, 9010, 9020, 9022, 9110, 9130, 9132 cardiofax Q, 9620 ...

