


## EE103BK-ML Leadwires for GE Marquette Multi-Link systems



NB the product might differ from the illustrations and it is subject to changes without prior notice.

### Description

Leadwires for GE Marquette Multi-Link.

The plastic connector towards the device is conformal to the original. The terminations have 4mm Banana Plugs, IEC coded.

Every trademark is propriety of the respective owner.

### Technical features

<b>Length:</b>	<ul style="list-style-type: none"> <li>✓ Chest leads: 73 cm</li> <li>✓ Limbs leads: 103 cm</li> </ul>
<b>Terminations:</b>	4mm Banana
<b>Yoke:</b>	not present
<b>Connector:</b>	bipolar, GE AAMI
<b>Color code:</b>	IEC
<b>Resistors:</b>	not present
<b>Protection against defibrillation:</b>	not present
<b>Packaging:</b>	No. 1 set



20%



80%

5° C  
41° F



40° C  
104° F

