


EE103BK-TC Leadwires for Philips electrocardiographs TC series



NB the product might differ from the illustrations and it is subject to changes without prior notice.

Description

Leadwires for PHILIPS electrocardiographs TC series.

The plastic connector towards the device is conformal to the original. The terminations have 4mm Banana Plugs, IEC coded.

Every trademark is propriety of the respective owner.

Technical features

Length:	<ul style="list-style-type: none"> ✓ Chest leads: 68 cm ✓ Limbs leads: 112 cm
Terminations:	4mm Banana
Yoke:	not present
Connector:	bipolar and multipolai, Philips proprietary
Color code:	IEC
Resistors:	not present
Protection against defibrillation:	not present
Packaging:	No. 1 set (2 pcs chest multi-leads, 4 pcs. Individual limbs leads)

