


EE200-xDG 3- or 5-LEADS TRUNK CABLE FOR PHILIPS MONITOR EE200-xDGGI 3- or 5LEADS LEADWIRES FOR PHILIPS MONITOR



EE200-3DG



EE200-3DGGI

NB the product might differ from the illustrations and it is subject to changes without prior notice.

Description

Trunk cable and leadwires for Philips monitor, series M1001B, M1002B, M3000A, M3001A, M3002A, MP2, M8102A, MP5, M8105A, 863063, 863064, 863065, 863066, 863068, 863077, 863085, 863086, 863087, 86308 etc. The terminations have Grabber connectors, IEC coded.

1 kOhm resistors are present, to record Breathing signals.

Every trademark is propriety of the respective owner.

Technical features

Leadwires	<ul style="list-style-type: none"> ✓ Connector: Philips 3- or 5-leads flat connector ✓ Terminations: Grabber
Trunk cable	<ul style="list-style-type: none"> ✓ Towards leadwires: Philips 3- or 5-leads flat connector ✓ Towards device: Philips round connector
Lenght	<ul style="list-style-type: none"> ✓ Leadwires 90 cm ✓ Trunk cable: 2,7 m
Packaging	No. 1 trunk cable or No. 1 set of leadwires
Compatibility	<p>Philips modules series M1001B, M1002B, M3000A, M3001A, M3002A, MP2, M8102A, MP5: M8105A.</p> <p>Philips monitor series SureSign 863063, 863064, 863065, 863066, 863068, 863077, 863085, 863086, 863087, 86308.</p>

