


EE203SK-AP6 CABLE FOR GE TELEMETRY APEX PRO & CH (6 leads)



NB the product is subject to changes without prior notice.

Description

6-leads cable for GE Telemetry CareScape PDM, GE Apex Pro and Apex Pro CH.
Snap terminations and palctic connector.
Unparalleled compatibility with the original cables. The connector makes a perfect match with the transmitter.

Every trademark is propriety of the respective owner.

Technical features

Length:	84 cm
Terminals:	Snap
Yoke:	not present
Connector:	GE
Color coding:	IEC
Resistors:	not present
Defi protection:	not present
Packaging:	1 cable each

