

EE203SK-IV3 PATIENT CABLE FOR PHILIPS INTELLIVUE 3 LEADS
EE203SK-IV5 PATIENTCABLE FOR PHILIPS INTELLIVUE 5 LEADS



NB the product might differ from the illustrations and it is subject to changes without prior notice.

Description

Patient cable for Philips telemetry Intellivue, 5- and 3-leads.
The plastic connector towards the device is conformal to the original, with the shower protection. The terminations have Snap connectors, IEC coded.

Every trademark is propriety of the respective owner.

Technical features

Lunghezza:	88 cm
Terminali:	Snap
Partitore:	not present
Connettore:	Philips DG 3- or 5-leads, with shower protection
Codifica colori:	IEC
Resistenze:	not present
Protezione da defibrillazione:	not present
Confezionamento:	No. 1 cable each

