



"EUROSPIRO" Disposable Blue Large Nose Clips



DESCRIPTION

"EUROSPIRO" Disposable Nose clips are intended for spirometry tests in order to reduce the risk of introducing a technical error. LUMED rubber- and foam-padded nose clips help prevent patients from inhaling or exhaling through the nose during spirometry.

LUMED nose clips are made of white thermoplastic, latex-free. The tension in the nose clip body keeps it firmly in place while the rubber and foam pads are designed for the best patient comfort.

The nose clips are disposable and intended for single-patient use only.

TECHNICAL FEATURES

REFERENCE	TSTRA070/10	TSTRA070/100	TSTRA090/10	TSTRA090/100
DESCRIPTION	Disposable Rubber- and Foam-padded Nose Clips			
NOSE CLIP BODY MATERIAL	Polyacetate			
NOSE CLIP PADS MATERIAL	Polyvinyl	Polyvinyl	Ethylene-vinyl acetate	Ethylene-vinyl acetate
COLOUR	Blue			
DIMENSIONS	70 x 45 mm			
MANUFACTURER	LUMED srl Via Staffora, 18/9 - 20073 Opera (MI) - ITALY			
SALE PACKAGING	10 pcs in cardboard box	100 pcs in cardboard box	10 pcs in cardboard box	100 pcs in cardboard box
STORAGE	Temperature: 0°C +50°C - Relative humidity: 20% + 80%			