


**EE100CRI-200 PATIENT CABLE FOR RESTING ECG**  
**EE100CRSI-200 PATIENT CABLE FOR STRESS-TESTING**  
**EE100CRSI-350 "XL" PATIENT CABLE FOR STRESS-TESTING**



NB the product is subject to changes without prior notice.

**Description**

One-piece patient cable for Stress-Testing ECG (12-leads). It is completely compatible with the original ones. The plastic connector towards the device is conformal to the original. The terminations have Snap plugs. Made for devices from the following Brands: Bosch, Cardioline, Dimeq, Kenz. The Compatibility List could be extended to further devices, according to the specific features.

Every trademark is propriety of the respective owner.

**Technical Features**

<b>Length:</b>	<ul style="list-style-type: none"> <li>✓ Trunk: <ul style="list-style-type: none"> <li>✓ EE100CRI-200, EE100CRSI-200: 200 cm</li> <li>✓ EE100CRSI-350: 350 mm</li> </ul> </li> <li>✓ Leadwires: chest 75 cm, limbs 100 cm</li> </ul>
<b>Terminations:</b>	<ul style="list-style-type: none"> <li>✓ EE100CRI-200: 4 mm Banana</li> <li>✓ EE100CRSI-200, EE100CRSI-350: Snap</li> </ul>
<b>Yoke:</b>	compact, sealed
<b>Connector:</b>	Kenz PC-104 type round 16-pin connector, made in plastic.
<b>Colour code:</b>	IEC
<b>Resistors:</b>	10 kΩ with defi protection
<b>Compatibility:</b>	BOSCH EKG 601. CARDIOLINE Delta 1, 3, 30, 30D, 60, 60D, 1plus, 3plus, 3/6plus. Elan, Elan1100. Delta 60plus, ElanStress. KENZ ECG-103, 104, 106, 107 ...

