


EE100ERI-200 PATIENT CABLE FOR RESTING ECG EE100ERSI-200 PATIENT CABLE FOR STRESS-TESTING



NB the product is subject to changes without prior notice.

Description

One-piece patient cable for Rest-ECG (12-leads diagnostic ECG). It is completely compatible with the original ones. The plastic connector towards the device is conformal to the original. The terminations have 4mm Banana plugs. The Compatibility List could be extended to further devices, according to the specific features.

Every trademark is propriety of the respective owner.

Technical features

Length:	✓ Trunk: 200 cm ✓ Leadwires: chest 70 cm, limbs 110 cm
Terminations:	✓ EE100NRI-200: 4mm Banana ✓ EE100NRSI-200: Snap
Yoke:	molded and sealed
Connector:	D-sub 15 pin, NK style. Made in plastic.
Colour code:	IEC
Resistors:	10 kΩ
Defi protection:	present
Compatibility:	BENEWARE CS-18x, -28x. BIO CARE ECG101, 101G, 300, 300G, 6010, 8080, 3030, 9801, 9803, 9803G, 1200, 3000, 3010, 3030, 6010, 1210, 1216, 1230 (12-channel), 9803. BIOLIGHT BLT-1203, BLT E12A, BLT E12, E70, E80. CAREWELL ECG-1101B/G, -1103B/G/W/L/LW, -1106L, -1112L/M, T10, T12. PCECG-500. CMICS-DONGJANG ECG-11A, 11B, 11D, 11D(II), 120, 1206, 1212, 1220, 1230S, 1230T, 32A, 32A(II) & all models with D-sub 15 connector. COMEN CM100, 300, 600, 1200, 1200A, 1200B. CONTEC all stand-alone models with D-sub15 connector, Contec 8000. COSMED Quark C12X. DMS DM SOFTWARE Cardioscan Resting 12-lead. EDAN SE-1, -100, -3, -3WS, -300, -301, -6, -600, -601A-B-C, -12, -12 express, -1200, -1200 express, -1010 (not wireless), -1201, -1515 DP12, ST-1212; VE-100, -300, -1010. EDAN-iNS ECG-1A, -1B, -3A, -3A plus, -3B, -3B plus, -6A, -6B, -12A. FAZZINI all models with D-sub15 connector. I-MEDIK EKG-312 touch, EKG-6012 touch. IMED Cardiax pc-ecg. INTERMED ECG-1101B/G, -1103B/G/W/L/LW, -1106L, -1112L/M, T10, T12. PCECG-500. KALAMED KES-301, -601, -121T, KES-1000. LABORATORIOS MACRIMASA Vitalfic E121, Vitalfic E123. M&B ECG-1203, -1206, -

1212. MINDRAY DECG-03A, BeneHeart R3, BeneHeart R12. MORETTI DIMED all models with D-sub15 connector, pc-ecg SEIVA. NORAV 1200-M, 1200-HR, 1200-HR-T, NHH-1200, 1200-S. NECG3+, NECG-12+, NECG-12C, NECG-12C+. NECG-3, NECG-12, NECG-12+. OSEN ECG-8130A, -8130, 8110, -8112. ECG Workstation OS700. PROGETTI all models with D-sub15 connector & MyECG, PG ECG. SINOHERO SE101, SE505, SE508, SE600, SE1200, SE1200 Light. SONOSCAPE all models with D-sub15 connector. SONOSTAR SE-1, -1A, -3, -3B, -3C, -6, -6A, 12A, -12B, -12P, Ebox-12. SPACELABS CardioExpress SL3, SL6, SL12. VALES & HILLS CV200. CV1200. ZONCARE ZQ-1201, -1201G, -1203, -1203G, -1203S, -1206, -1212, IMAC730, T12 ...

