


EE100SRI-200 patient cable for Resting ECG
EE100SRSI-200 patient cable for Stress-Testing
EE100SRSI-350 "XL" patient cable for Stress-Testing, Snap
EE100SRGI-350 "XL" patient cable for Stress-Testing, Grabber



NB the product is subject to changes without prior notice.

Description

One-piece patient cable with 12 leadwires, for Resting and Stress-ECG. Cables from EE100SRxl-xx0 series are 100% compatible. The connector towards device is made in plastic and it is conformal to the original. The terminations are coded with IEC colors scheme and have these connectors: 4mm Banana (EE100SRI-200), Snap (EE100SRSI-200 and -350) or Grabber (EE100SRGI-350). The Compatibility List could extend to further devices, according to the specific features.

Every trademark is propriety of the respective owner.

Technical data

Length:	Main cable: ✓ EE100SRI-200, EE100SRSI-200: 200 cm ✓ EE100SRSI-350, EE100SRGI-350: 350 mm Leadwires: chest 75 cm, limbs 100 cm
Terminations:	✓ EE100SRI-200: 10 leads 4mm Banana ✓ EE100SRSI-200, EE100SRSI-350: 10 leads Snap ✓ EE100SRGI-350: 10 leads Grabber
Yoke:	molded, compact style
Connector:	Cannon style, standard 15 pin, plastic casing with screws
Colour code:	IEC
Protection:	present, with 10 kΩ resistors
Compatibility:	ESAOTE P80, P80six, P80power, P8000, P8000power, P8000 tele, Bluecg, MyCardiopad, MyCardiopad XL, AT-10 plus, AT-102 plus, Cardiovit M2015 etc. FUKUDA DENSHI P80. KONTRON EKG P80. SCHILLER AT-1, AT-2, AT-2plus, AT-2 light, AT-4, AT-101, AT-104, AT-110. CS-6, CS-6/12, CS-12, AT-5, AT-6, AT-10, AT-10+, AT-60, AT-100. MS-2007, MS-2010, MS-2015. AT-10plus, AT-102, AT-104PC, CS-100, CS-200, CS200 excellence, CS200 touch etc. SECA Serie CT3000i, CT6i/6Pi, CT8000i/P/L/T. CT110, CT220. CT cardiopad; PT-160. WELCH ALLYN Serie AT. Serie SP. Serie CT, CP50, CP150, CP300. ZIMMER Cardio 5, Cardiacompact ...

