



## EE103BK-CAM Leadwires for GE Marquette CAM 14 and HD



NB the product might differ from the illustrations and it is subject to changes without prior notice.

### Description

Leadwires for ECG modules GE Marquette CAM14 and HD.  
The plastic connector towards the device is conformal to the original. The terminations have 4mm Banana Plugs, IEC coded.

Every trademark is propriety of the respective owner.

### Technical features

<b>Length:</b>	<ul style="list-style-type: none"> <li>✓ Chest leads: 73 cm</li> <li>✓ Limbs leads: 103 cm</li> </ul>
<b>Terminations:</b>	4mm Banana
<b>Yoke:</b>	not present
<b>Connector:</b>	bipolar, GE AAMI
<b>Color code:</b>	IEC
<b>Resistors:</b>	not present
<b>Protection against defibrillation:</b>	not present
<b>Packaging:</b>	No. 1 set of 10 leadwires.

