


## EE103BK-TC16 set of 16 Leadwires for Philips TC series



NB the product might differ from the illustrations and it is subject to changes without prior notice.

### Description

Leadwires for PHILIPS electrocardiographs TC series. Set of 16 leads: N, R, L, F, C1-C6, V3R, V4R, V5R, V7, V8, V9

The plastic connector towards the device is conformal to the original. The terminations have 4mm Banana Plugs, IEC coded.

Every trademark is propriety of the respective owner.

### Technical features

<b>Length:</b>	<ul style="list-style-type: none"> <li>✓ Chest leads: 68 cm</li> <li>✓ Limbs &amp; additional leads: 112 cm</li> </ul>
<b>Terminations:</b>	4mm Banana
<b>Yoke:</b>	not present
<b>Connector:</b>	bipolar and multipolar, Philips proprietary
<b>Color code:</b>	IEC
<b>Resistors:</b>	not present
<b>Defibrillation protection:</b>	not present
<b>Packaging:</b>	2 pcs chest multi-leads, 10 pcs. Individual leads (4 pcs chest, 3 pcs right chest, 3 pcs posterior)

