



"EUROCATCH" ECG Suction Electrodes



DESCRIPTION

LUMED "EUROCATCH" Reusable Precordial ECG Suction Chest Electrodes have features that make them reliable during signal recording with all ECG devices available on the market. The metallic part is silver and silver chloride (ag/agCl) in order to allow a safe and clear signal transmission.

LUMED reusable precordial ECG suction chest electrodes, thanks to screw and clip connections, can be used with all ECG devices and monitors with patient cable that has banana pin (4mm) or snap. The product is latex and nickel free. They are available in set of 6 suction electrodes, as showed in the picture, and in the following versions: adult (EL0002) and paediatric (EL0001).

TECHNICAL FEATURES

REFERENCE	EL0001 <i>paediatric</i>	EL0002 <i>adult</i>
DESCRIPTION	Reusable precordial ECG suction chest electrodes	
ELECTRODE MATERIAL-COATING	SILVER/SILVER CHLORIDE (ag/agCl)	
CONNECTION	UNIVERSAL - hole with screw (4 mm) and clip	
BALL MATERIAL	PLASTIC (MEDICAL GRADE PVC)	
METALLIC PART DIAMETER	15 mm	24 mm
BALL DIAMETER	25 mm	33 mm
COMPLIANCE TO STANDARDS	Directive 93/42/EEC and amendments. Class I medical device Italian Medical Device Data standard Registration: 1576725; CND: C020501; CE: 0476 *Non-Active Medical Device*	
MANUFACTURER	LUMED srl Via Staffora, 18/9 - 20090 Opera (MI) - ITALY	
SALE PACKAGING	6 pcs in polyethylene pouch	
STORAGE	Temperature: 0°C + +50°C - Relative humidity: 20% + 80%	
PRODUCT SHELF LIFE	5 years from production date printed on packaging	