

BMSY013 – THERMAL PAPER FOR SONY VIDEOPRINTERS



Description

LUMED SRL offers a new model of high-quality compatible film for Sony videoprinters (**UPP110-HD**). Thanks to the optimal grey-scale response, our product solved the major problem of the other compatible films, because it can catch all the clinical information, like as the original does.

End-users of BMSY013 also benefit from relevant saving, due to the lower price and the high quality of our film.

Each roll is packed individually, with Aluminum foil. Each box contains 25 rolls.

Technical features

REF:	BMSY013
DESCRIPTION:	thermal fil for videoprinter
SIZE:	roll height: 110 mm – roll length: 20 m
PACKAGING:	Aluminum film
STORAGE:	Temperature: < 40°C – Relative Humidity: < 65%
	Do not expose the film to reactive gases, sun light, UV sources etc.

Directions for use

BMSY013 has the same directions for use of the original film is. Follow the directions for storage available on the package of each roll. Do not expose the film to reactive gases, sun light, UV sources etc. Moreover, avoid any contact with carbon paper, plastifying agents, self-adhesive ribbons, alcohol, Ketons, Esters and organic substances containing Esters, cosmetics, and products for skin care. Do not touch BMSY013 with wet or moistened hands.

Conformity

Directive 93/42/CEE Class I MDD