



CE



# EUROCATCH

IS MADE IN LUMED®








**LUMED®**  
www.lumed.com



## EUROCATCH È UN PRODOTTO LUMED® / IS MADE IN LUMED®

Elettrodi monouso a gel solido o liquido per ECG, prova da sforzo, Holter e monitoraggio.

Solid & liquid Gel disposable electrodes, for Resting ECG, Stress-Test, Holter and Monitoring.

										
REF	ET1166	ET4166	ET1119	ET1136E	ET1184	ET1140	ET4140			
USO PURPOSE	ECG a riposo Resting ECG		Holter Monitor	Stress-Test						
DIMENSIONI SIZE	Quadrato / Square 32 x 32 mm		Rotondo / Round 50 mm	Ovale / Oval 36x50 mm						
GEL	Liquido Liquid	Solido Solid	Liquido Liquid	Liquido Liquid	Liquido Liquid	Liquido Liquid	Liquido Liquid	Solido Solid		
SUPPORTO SUPPORT	Spugna Foam	Spugna Foam	Spugna Foam	Spugna Foam	Spugna Foam	Spugna Foam	Spugna Foam	Spugna Foam		
ALTRO OTHER	-	-	-	-	Connettore Connector Banana 4 mm	Radiotrasparente, amagnetico Radiotransparent, amagnetic				

Tutti i marchi registrati o depositati sono proprietà dei rispettivi titolari. I prodotti potrebbero differire da quanto illustrato. Il produttore si riserva il diritto di modificare i prodotti qui descritti senza darne preavviso.

Every Trademark or Registered Mark is property of the respective owner. The products may differ from the content of the present document. The Manufacturer has the right to modify the products without prior notice.



LUMEDsrl  
info@lumed.com

VIA STAFFORA 18/9 • 20090, OPERA [MI] • ITALY  
TEL. +39 02 57606751 • FAX +39 02 57606752

VIA SENIO 36/40 • 47121, FORLÌ [FC] • ITALY  
TEL. +39 0543 702380 • FAX +39 0543 707518