



"EUROCARD" Universal Medical Cart



N.B. Product is subject to change without prior notice. The actual product might differ from the pictures here shown. Arm and additional shelf are optional accessories.

DESCRIPTION

LUMED "EUROCARD" Universal Medical Cart is a product line for equipment intended for medical ambulatory and clinics. 1.5 mm thick steel structure that guarantees durability and reliability for customers. Several extra accessories are available. Especially, the articulated arm does not have a simple cable holder function but it allows a complete coverage thanks to the reclosable snap lids with an easy closure. On the rear panel there is a compartment for daily use accessories (paper, electrodes, gel) and to support patient cable leads with reusable electrodes while you are using it. It is painted in order to be scratch resistant and very easy to clean. It is equipped with four 70 mm diameter wheels and, two of them, with brake. The EUROCARD is tested for stability on an inclined plane up to 12°. It is supplied unassembled (assembly kit includes the necessary tools and instructions).

TECHNICAL FEATURES

REFERENCE	EP-LU4003	EP-LU4006
DESCRIPTION	Trolley, simple version with base and 4 wheels (2 with brakes), shelf and rear panel with compartment.	
COLOUR	Single color light-grey RAL 7035	light-grey RAL 7035 and blue RAL 5012 (shelf)
DIMENSIONS	40 x 45 cm (L x W) 90 cm (total height from pavement)	
WEIGHT	around 15 kg	
COMPLIANCE TO STANDARDS	Directive 93/42/EEC and amendments. Class I medical device Italian Medical Device Data standard Registration: 1241013/R; CND: V 9099; GMDN: 34909; CE: 0476	
MANUFACTURER	LUMED srl Via Staffora, 18/9 - 20090 Opera (MI) - ITALY	
PAYLOAD	15 kg each shelf; 40 kg max capacity	
EXTRA ACCESSORIES	<ul style="list-style-type: none"> - EP-LU4002 arm for cable (RAL 7035) - EP-LU4004 additional shelf (RAL 7035) - EP-LU4005 storage unit with lid (RAL 7035) - EP-LU4017 support for US probes (RAL 7035) 	<ul style="list-style-type: none"> - EP-LU4007 arm for cable (RAL 7035, lids RAL 5012) - EP-LU4008 additional shelf (RAL 5012) - EP-LU4009 storage unit with lid (RAL 7035, lid RAL 5012) - EP-LU4017 support for US probes (RAL 7035)