

EP-LU6000 "EUROCARD" trolley for medical systems



N.B. Product is subject to change without prior notice. The actual product might differ from the pictures here shown. Arms, monitors, computer and printer are optional accessories..

DESCRIPTION

LUMED EUROcard 6000 Medical Trolley for PC-based systems, destined for ambulatories and clinics. 1.5 mm thick steel structure that guarantees durability and reliability for customers. Several extra accessories are available. 4 swivel castors. Opening on the front and rear door with lock. Hinges can be mounted on left or right side. 4 shelves\surfaces available. EUROcard 6000 is powder-painted in order to be scratch resistant and very easy to clean.

FEATURES

DESCRIPTION	Base with 4 swivel castors (2 with brake, 2 antistatic), body, top, rear door with lock, 2 internal shelves.
COLOUR	Light-grey RAL 7035 and blue RAL 5012
SIZE	<ul style="list-style-type: none"> - Body 37,5 x 65 cm (D x W) - Base 62 x 72 cm (D x W) - Total height approx. 105 cm
WEIGHT	Approx. 30 Kg
LOAD	35 Kg (single shelf); 70 Kg (overall load capacity)
ACCESSORIES	<ul style="list-style-type: none"> - Keyboard support - EUROvacu vacuum °ECG patient cable for Resting ECG and Stress-Test - Hinged arm for patient cable - Hinged arm for probe - Rack for probes - Support for 3 liters calibration syringe - Arm for monitor (with VESA support) - <i>EUROcard6000 can mount up to 3 arms for monitor</i> - Basket - Isolation transformer