




"EUROSPIRO" Disposable Paper Mouthpieces



DESCRIPTION

LUMED "EUROSPIRO" Disposable Cardboard Mouthpieces are made from durable Kraft paper and designed for single patient use. Each mouthpiece is individually wrapped in order to guarantee the maximum hygiene levels. They are available in a variety of sizes to fit most pulmonary function equipment and are packaged in convenient dispenser boxes of 100 or 150 pieces or bulk cases of 500 pieces. LUMED "EUROSPIRO" Disposable Cardboard Mouthpieces meet the requirements of the applicable EC directives in order to obtain the CE marking.

TECHNICAL FEATURES

REFERENCE	SPIROMETER	M/F	INTERNAL DIAMETER	EXTERNAL DIAMETER	LENGTH	QUANTITY
TS0048	JAEGER	M	29,7 ±0,1	32,0 ±0,5	70 ±2	standard: 500 dispenser: 4 x 100*
TS3456	Sibelmed	M	22,2 ±0	25,0 -0; +0,2	70 ±2	standard: 500 dispenser: 4 x 100*
TS3962	MIR, Vitalograph, Micromedical	M	28,0 ±0,2	30,0 -0; +0,5	70 ±2	standard: 500 dispenser: 4 x 100*
TS3963	Fukuda Sangyo, Medical Graphics, Biomedin, Custo-Vit	M	25,3 ±0,2	27,0 -0; +0,5	70 ±2	standard: 500 dispenser: 4 x 100*
TS5218	COSMED	M	28,0 ±0,2	30,5 ±0,2	70 ±2	standard: 500 dispenser: 4 x 100*
TS5422	Custo-Vit M	F	29,2 +0; -0,2	32,0 -0; +0,5	70 ±2	standard: 500 dispenser: 4 x 100*
TS8182	SensorMedics, Microgard	F/(M)	30,3 ±0,2	33,0 -0; +0,5	70 ±2	standard: 500 dispenser: 4 x 100*
TS8248	Sibelmed, BTL	M	26,2 ±0,2	28,0 -0; +0,5	70 ±2	standard: 500 dispenser: 4 x 100*
TS8551	Universal Pediatric Size	F/M	18,5 ±0,2	22,5 -0; +0,5	70 ±2	standard: 500 dispenser: 4 x 150*
TS8822	VB filters	M	25,4 ±0,2	26,8 -0; +0,5	70 ±2	standard: 500 dispenser: 4 x 100*
TS9190	Sorin LifeWatch, Fukuda Denshi, Welch Allyn	M	24,2 ±0,2	26,0 -0; +0,5	70 ±2	standard: 500 dispenser: 4 x 100*

*TSxxxx/1

Note: To order the cardboard mouthpieces in dispenser, please specify "/1" at the end of the reference

MOUTHPIECE MATERIAL	Kraft Paper
WRAPPING MATERIAL	Neutral Polypropylene Film
COMPLIANCE TO STANDARDS	Directive 93/42/EEC and amendments. Class IIa medical device Italian Medical Device Data standard Registration: 49378/R; CND: R 9001; GMDN: 44545; CE: 0476 *Non-Active Medical Device*
MANUFACTURER	LUMED srl Via Staffora, 18/9 - 20090 Opera (MI) - ITALY
STORAGE	Store in a cool and dry place