

DESCRIZIONE / DESCRIPTION	TIPO / TYPE	TERMINAZIONI / TERMINATIONS		
		BANANA 4 MM	SNAP	SNAP XL (L = 4.26 m)
AMBISEA AV-9100, 9300, 9600, 9120.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
AMEDTEC ECGpro, CardioPart 12 Blue, CardioPart 12 USB.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
ASPEL AsCard 3 (1992 – 1993).	completo, one-piece	EE100CUBI-200	EE100CUSI-200	
ASTRON 3/12b, 3/12i; Comp1x12, Comp3x12, Comp6x12, Comp3, Comp12, Comp60.	completo, one-piece	EE100DRI-200	EE100DRSI-200	EE100DRSI-350
ASTRON Replay, Replay 3.	completo, one-piece	EE100NRI-NV200		
ATES EBNEURO Easy ECG Pocket & Mobile.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
AVECINNA Sina 100.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
AVIONICS DEL MAR 910, 920, 930, 940, 942, 950, 960, 970, 980.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
BENEWARE CS-18x, -28x.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
BIOCARE ECG101, 101G, 300, 300G, 6010, 8080, 3030, 9801, 9803, 9803G, 1200, 3000, 3010, 3030, 6010, 1210, 1216, 1230 (12-channel), 9803.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
BIOLIGHT BLT-1203, BLT E12A, BLT E12, E70, E80. CAREWELL ECG-1101B/G, -1103B/G/W/L/LW, -1106L, -1112L/M, T10, T12. PCECG-500.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
BIONET cardiocare 2000, cardiotech 3000, cardio 7.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
BMS BIOMEDICAL SYST. 5 DIN-Snap.	terminazioni, leadwires		EEHL05-60	
BMS BIOMEDICAL SYSTEMS 7 DIN-Snap.	terminazioni, leadwires		EEHL07-60	
BOSCH 501A, 501NV, 503, 5038, 503D, 503NV, 503V, 506, 603D, 603ER, 606.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
BOSCH 501A, 501NV, 503, 503B, 503D, 503NV, 503V, 506, 506D, 601.	completo, one-piece	EE100NRI-NV200		
BOSCH EKG 601	completo, one-piece	EE10CRI-200	EE10CRSI-200	EE100CRSI-350
BOSCH 603D, 603ER, 603S, 606.	completo, one-piece	EE100CUBI-200	EE100CUSI-200	
BRAEMAR 5x DIN-Snap.	terminazioni, leadwires		EEHL05-60	
BRAEMAR DL-400, DL-450, DL-700, DXP-1000.	terminazioni, leadwires		EEHL07-60	
BTL Serie BTL-08; CardioPoint-ECG C300, C600. CardioPoint- ERGO C300, C600.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350

DESCRIZIONE / DESCRIPTION	TIPO / TYPE	TERMINAZIONI / TERMINATIONS		
		BANANA 4 MM	SNAP	SNAP XL (L = 4.26 m)
BURDICK EK10, Elite, E350, E350I, E550, E560, Eclipse LE, LE I, LE II, Elite II.	completo, one-piece	EE100NRI-NV200		
BURDICK QUINTON Eclipse 800, Eclipse LE/LEi (> SN 1051), CareCenter MD.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
CAMBRIDGE VS-6, SC1000.	completo, one-piece	EE100NRI-NV200		
CAMBRIDGE CM6000, CM6400.	completo, one-piece	EE100CUBI-200	EE100CUSI-200	
CAMEDIO 12B, 12H, 12S.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
CAMINA CARDIORAPID K11, K111, K112pocket, K130B, K130D, K131, K30P, K300P, K31, K312P, K33M, K36M, K360B, K360D, K3600LCD, K600P, K600M, K612.	completo, one-piece	EE100DRI-200	EE100DRSI-200	EE100DRSI-350
CAMINA CARDIOSWIFT K112LCD, K36LCD, conn. D-sub 15.	completo, one-piece	EE100DRI-200	EE100DRSI-200	EE100DRSI-350
CARDIAC SCIENCE X12+.	completo, one-piece		EEHL20XL	
CARDIETTE, ELETTRONICA TRENTINA Hr One, Hr One48 5x DIN.	terminazioni, leadwires		EEHL05-60	
CARDIETTE - ELETTRONICA TRENTINA HR one 48 7x DIN, DXP1000, Microscan 2000.	terminazioni, leadwires		EEHL07-60	
CARDIETTE - ELETTRONICA TRENTINA every version of AR600,1200, 2100; Start 100P, Start 100H, Start 200, 200HV. DOC 101, 103; Autoruler 12/1, 12/3; Excel 103, 106, 106SCP; Excel 112, 212; Daedalus 346; Serie Daedalus View.	completo, one-piece	EE100DRI-200	EE100DRSI-200	EE100DRSI-350
CARDIETTE - ELETTRONICA TRENTINA every model of Micro, Microruler, Microtel.	completo, one-piece	EE100NRI-NV200		
CARDIO CONTROL Cardio Perfect.	completo, one-piece	EE100CUBI-200	EE100CUSI-200	
CARDIOLEX EC Sense.	completo, one-piece	EE100CUBI-200	EE100CUSI-200	
CARDIOLINE all versions of AR600,1200, 2100; E1, E1R, E3; ECG100+, ECG200+, ECG100S, ECG200S.	completo, one-piece	EE100DRI-200	EE100DRSI-200	EE100DRSI-350
CARDIOLINE Click Holter.	terminazioni, leadwires		EEHL05-60	
CARDIOLINE Delta 1, 3, 30, 30D, 60, 60D, 1plus, 3plus, 3/6plus. Elan, Elan1100. Delta 60plus, ElanStress	completo, one-piece	EE10CRI-200	EE100CRSI-200	EE100CRSI-350
CARDIOLINE Delta 300, 360.	completo, one-piece	EE100CUBI-200	EE100CUSI-200	
CARDIOLINE E1, E1R, E3, Delta300, Delta360, Delta1, Delta3.	completo, one-piece	EE100NRI-NV200		
CARDIOMQA 100H, 300H, 600H, 1200H.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
CARDIONICS Cardio Touch, Cardio Plug.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350

DESCRIZIONE / DESCRIPTION	TIPO / TYPE	TERMINAZIONI / TERMINATIONS		
		BANANA 4 MM	SNAP	SNAP XL (L = 4.26 m)
CARDIONICS conn. D-sub15.	completo, one-piece	EE100NRI-NV200		
CARDIOTEK serie EP-Tracer.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
CAREFUSION JAEGER Masterscreen ECG.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
CAREWELL CHS-12 Dynamic ECG.	terminazioni, leadwires		EEHL10-60	
CMICS–DONGJANG ECG-11A, 11B, 11D, 11D(II), 120, 1206, 1212, 1220, 1230S, 1230T, 32A, 32A(II) & conn. D-sub 15.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
COLIN YM406i.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
COLSON DUPONT Cardi-3, E1, E3. Cardiopratique. Models from these series: AR, Cardipocket, Cardio, Excel, Autoruler, Microruler; Daedalus 346.	completo, one-piece	EE100DRI-200	EE100DRSI-200	EE100DRSI-350
COLSON DUPONT every model of Microruler.	completo, one-piece	EE100NRI-NV200		
COMEN CM100, 300, 600, 1200, 1200A, 1200B.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
CONTEC stand-alone models conn. D-Sub 15, Contec 8000.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
CONTEC TLC4000, TLC5000, TLC6000.	terminazioni, leadwires		EEHL10-60	
CONTEC TLC49804.	terminazioni, leadwires		EEHL07-60	
CONTEC TLC9803, PM85, PM60B.	terminazioni, leadwires		EEHL05-60	
CORSCIENCE BT12 wireless.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
COSMED Cardico 306, 601H, 1201, 1211.	completo, one-piece	EE100KRI-200		
COSMED Quark C12, T12.	completo, one-piece		EEHL20XL	
COSMED Quark C12X.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
CUSTO.MED card M, cardio 130.	completo, one-piece	EE100CUBI-200	EE100CUSI-200	
CUSTO.MED norm, cardio130.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
DAHIAN Yasham 101, 310; Araz 630.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
DATASCOPE Spectrum, Passport 2.	completo, one-piece		EEHL20XL	
DATASCOPE DT-4500	completo, one-piece		EE203SK-APFH	

DESCRIZIONE / DESCRIPTION	TIPO / TYPE	TERMINAZIONI / TERMINATIONS		
		BANANA 4 MM	SNAP	SNAP XL (L = 4.26 m)
DATRIX INTRICON 7 DIN-Snap.	terminazioni, leadwires		EEHL07-60	
DATRIX INTRICON XR-300 and others with 5x DIN-Snap.	terminazioni, leadwires		EEHL05-60	
DAVINSA models conn. D-Sub 15	completo, one-piece	EE100ERI-200	EE100ERSI-200	
DEGO Cardiomed 3.	completo, one-piece	EE100CUBI-200	EE100CUSI-200	
DEL MAR Digicorder 419, 420, 445, 447, 448, 449, 453, 456 (red, white, green, black, brown).	terminazioni, leadwires		EEHL05-60	
DEL MAR Digicorder 423, 423A, 459, 461, 461A, 463, 464, 483, 485 (white, black, red, green, orange, blue, brown).	terminazioni, leadwires		EEHL07-60	
DEL MAR REYNOLDS SPACELABS Aria, Evo.	terminazioni, leadwires		EEHL40XL	
DMS DM SOFTWARE Cardioscan Resting 12-lead.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
DMS DM-SOFTWARE, MEIGAOYI 300-5, 300-540, 300-6, 300- 7, 300-7P, 300-7L, 300-7M.	terminazioni, leadwires		EEHL05-60	
DMS DM SOFTWARE - MEIGAOYI 300-540, 300-6, 300-7, 300-7P, 300-7L, 300-7M.	terminazioni, leadwires		EEHL07-60	
DMS DM-SOFTWARE, MEIGAOYI 300-8, 300-12L.	terminazioni, leadwires		EEHL10-60	
DR. LEE 310A, 310B.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
DR. LEE 120A, 310A, 310B.	completo, one-piece	EE100NRI-NV200		
DR. SCHMIDT SmartScript.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
DR. VETTER pc-ekg.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
EDAN SE-1, -100, -3, -3WS, -300, -301, -6, -600, -601A-B-C, - 12, -12 express, -1200, -1200 express, -1010 (not wireless), - 1201, -1515 DP12, ST-1212; VE-100, -300, -1010.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
EDAN-iNS ECG-1A, -1B, -3A, -3A plus, -3B, -3B plus, -6A, -6B, - 12A.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
ERGOLINE Autoruler; Autoscript 1/12, 3/12, 6/12; Ergoscript EK301Z, EK1012, EK3003, EK3006, EK3012.	completo, one-piece	EE100DRI-200	EE100DRSI-200	EE100DRSI-350
ESAOTE Aria, Evo	terminazioni, leadwires		EEHL40XL	
ESAOTE Lifecard CF 3-ch, 3 canali	completo, one-piece		EEHL30	
ESAOTE Lifecard CF 3-ch w/ belt clip, 3 canali con clip cintura	completo, one-piece		EEHL30/CLIP	
ESAOTE Lifecard CF 3-ch w/7x DIN, 3 canali 7x DIN	terminazioni, leadwires		EEHL07-60	

DESCRIZIONE / DESCRIPTION	TIPO / TYPE	TERMINAZIONI / TERMINATIONS		
		BANANA 4 MM	SNAP	SNAP XL (L = 4.26 m)
ESAOTE P80, P80six, P80power, P8000, P8000power, P8000 tele, Bluecg, MyCardiopad, MyCardiopad XL, AT-10 plus, AT-102 plus, Cardiovit M2015	completo, one-piece	EE100SRI-200	EE100SRSI-200	EE100SRSI-350
FAZZINI conn. D-Sub 15.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
FIDELITY LP3000.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
FUKUDA DENSHI Autocardiner: FCP11, FCP2101A, FCP220A, FCP4101A. FX101. MCP302. Cardimax: FX2111, FCP2155, FCP7101, FX3010, FX4010, FX7102, FX7202, FX7302, FX7402, FX7542, FX8222, FX8322, FX8328.	completo, one-piece	EE100FRI-200	EE100FRSI-200	
FUKUDA DENSHI Autocardiner: FCP11, FX101. MCP302. Cardimax: FX2111, FCP2155, FX4010, FX7402. Alpha 500, Alpha 600, FCP 2201.	completo, one-piece	EE100NRI-NV200		
FUKUDA DENSHI P80.	completo, one-piece	EE100SRI-200	EE100SRSI-200	EE100SRSI-350
FUKUDA ME Alpha 1000, 6600A, 6600AX.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
FUKUDA ME Cardiosunny 1000, 6000A, 6600AX.	completo, one-piece	EE100NRI-NV200		
GE MARQUETTE CAM14, CAM HD Set of 10 leadwires 4 mm Banana. Set 10 terminazioni banana 4 mm.	terminazioni, leadwires		EE103BK-CAM	
GE MARQUETTE Multi-Link Set of 10 leadwires 4 mm Banana. Set 10 terminazioni banana 4 mm.	terminazioni, leadwires		EE103BK-ML	
GE MARQUETTE Apex Pro FH, GE DT-4500.	completo, one-piece		EE203SK-APFH	
GE MARQUETTE Apex, Apex S	completo, one-piece		EE203SK-AP5	
GE MARQUETTE CareScape PDM, Apex Pro	completo, one-piece		EE203SK-AP6	
GE MARQUETTE 2-channel recorders serie 7000, 8000, Seer, Seer XT, Seer MC.	terminazioni, leadwires		EEHL05-60	
GE MARQUETTE HELLIGE 3-channel recorders serie 8000, Seer, Seer MC.	terminazioni, leadwires		EEHL07-60	
GE MARQUETTE HELLIGE Microsmart, Cardiosmart. MAC 400, 500, 600, 800, 1000, 1100, 1200, 1600, 2000, 3500	completo, one-piece	EE100GRI-200		
GEM-MED Serie Heart One +.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
GIMA Cardiopocket, Smart SE-1, -3, -6, -12 Express. Gimaplay 1. Cardiogima 1M, 3M, 6M, 12, 12M, Colour ECG. PC-ECG 33320, 33321; Acquisition box series MGY.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
GIMA ECG Holter ref. 35130.	terminazioni, leadwires		EEHL10-60	
HELLIGE EK36, EK41, EK43, EK56, EK403, EK413, EK512, EK512 E, EK512 P.	completo, one-piece	EE100CUBI-200	EE100CUSI-200	

DESCRIZIONE / DESCRIPTION	TIPO / TYPE	TERMINAZIONI / TERMINATIONS		
		BANANA 4 MM	SNAP	SNAP XL (L = 4.26 m)
HELLIGE EK51, Corina.	completo, one-piece	EE100NRI-NV200		
HELLIGE EK53.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
HEWLETT PACKARD PageWriter 100, 200, 300 (M1770A, M1771A, M1772A).	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
HONEYWELL PHILIPS Cardioplan ER330, CR300, CR330.	completo, one-piece	EE100NRI-NV200		
HORMANN VON BERG Bioset 3500, 3600, 3700, 6000, 8000.	completo, one-piece	EE100CUBI-200	EE100CUSI-200	
HORMANN VON BERG Bioset 9000.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
IMED Cardiax pc-ecg.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
I-MEDIK EKG-312 touch, EKG-6012 touch.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
INNOMED Serie Heart Mirror & Heart Screen. Cardio PC/E.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
INTERMED ECG-1101B/G, -1103B/G/W/L/LW, -1106L, -1112L/M, T10, T12. PCECG-500.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
KALAMED KES-301, -601, -121T, KES-1000.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
KENZ 108 & 110.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
KENZ 103, 104, 106, 107	completo, one-piece	EE10CRI-200	EE100CRSI-200	EE100CRSI-350
KENZ 108, 109, 110; Cardico 301, 302, 303, 305, 306, 601, 1201, 1207, 1207 STRESS, 1210, 1211, 1215	completo, one-piece	EE100KRI-200		
KONTRON EKG P80.	completo, one-piece	EE100SRI-200	EE100SRSI-200	EE100SRSI-350
LABTECH EC-12RT.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
M4MEDICAL M-trace, M-trace Mini, M-trace pc-ecg.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
MACRIMASA Vitalfic E121, Vitalfic E123. M&B ECG-1203, -1206, -1212.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
MAQUAIRE CT-100, -200, -200i, -200CV.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
MARQUETTE HELLIGE CardioSmart, CardioSys, Microsmart, MAC 500, MAC 1000, MAC 1200, Corina.	completo, one-piece	EE100CUBI-200	EE100CUSI-200	
MEDEA MedeaLine 3, MedeaLine 6p, Kardio PCM-u, Excercise Test, Cardio Test.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
MEDIANA YM303i, 406i, 412i, 812i.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350

DESCRIZIONE / DESCRIPTION	TIPO / TYPE	TERMINAZIONI / TERMINATIONS		
		BANANA 4 MM	SNAP	SNAP XL (L = 4.26 m)
MEDICAL ECONET Cardio-M, Cardio-M plus, Cardio-M PC, (instead of Quickels suction system / al posto del sistema a suzione Quickels).	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
MEDIGATE ECG-1211-DG; MeCA303i, 406i, 412i, 812i.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
MEDILOG OXFORD SCHILLER FD5.	terminazioni, leadwires		EEHL10-60	
MEDISPO EM9530, 9535, 9540, 8033, 9550, 9560.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
MEDSET Mobile Resting-ECG (Nihon-Kohden 1250, 1350, 1500, 1550).	completo, one-piece	EE100NRI-200	EE100NRSI-200	
MEDSET Serie Flashlight D-sub15 connector.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
MEGOS 3/12b, 3/12i; Comp1x12, Comp3x12, Comp6x12, Comp3, Comp12, Comp60.	completo, one-piece	EE100DRI-200	EE100DRSI-200	EE100DRSI-350
MEGOS Replay, Replay 3.	completo, one-piece	EE100NRI-NV200		
MESA Cardiax pc-ecg.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
MINDRAY Datascope Spectrum & Passport 2.	completo, one-piece		EEHL20XL	
MINDRAY DECG-03A, BeneHeart R3, BeneHeart R12.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
MORETTI DIMED models conn. D-Sub 15, pc-ecg SEIVA.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
MORTARA Eli 10, 10 mobile, 150, 210, 250, 350 (no WAM, no AM12).	completo, one-piece	EE100MRI-200	EE100MRSI-200	EE100MRSI-350
MORTARA AM12, WAM Set of 10 leadwires 4 mm Banana. Set 10 terminazioni banana 4 mm.	terminazioni, leadwires		EE103BK-WAM	
MORTARA Eli100, 200.	completo, one-piece	EE100NRI-NV200		
MORTARA X-12, H-12, X-12+, H-12+ (no T12, T12+); Surveyor II / X-12 Transmitter.	completo, one-piece		EEHL20XL	
MS-WESTFALIA LunoCard S & HR and models conn. D-Sub 15.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
NEUROSOFT Serie Poly-Spectrum: -12/E, -8/EX, -8/E, -8,-Velo.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
NIHON KHODEN ECG-1150 cardiofax C, 1250 cardiofax S, 1350 cardiofax M, 1500/-1550 cardiofax V, 2350 cardiofax M, PEA-1110. ECG-651, 9010, 9020, 9022, 9110, 9130, 9132 cardiofax Q, 9620.	completo, one-piece	EE100NRI-200	EE100NRSI-200	
NIHON KOHDEN Cardiofax 6151, 6353, 8110, 8240, 8270, 8330, 8340, 8350, 8370, 8420, 8423, 8450, 8453, 8820, 8823, 8830, 9020, 9022, 9022K, 9130, 9132. Cardiofax Q.	completo, one-piece	EE100NRI-NV200		

DESCRIZIONE / DESCRIPTION	TIPO / TYPE	TERMINAZIONI / TERMINATIONS		
		BANANA 4 MM	SNAP	SNAP XL (L = 4.26 m)
NORAV 1200-M, 1200-HR, 1200-HR-T, NHH-1200, 1200-S. NECG3+, NECG-12+, NECG-12C, NECG-12C+. NECG-3, NECG-12, NECG-12+.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
NORTH-EAST MONITORING 5 DIN-Snap.	terminazioni, leadwires		EEHL05-60	
NORTH-EAST MONITORING 7 DIN-Snap.	terminazioni, leadwires		EEHL07-60	
OSEN ECG-8130A, -8130, 8110, -8112. ECG Workstation OS700.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
OXFORD MEDILOG Excel, FD-4.	terminazioni, leadwires		EEHL07-60	
PHILIPS HP PageWriter 10, 10i, Pagewriter 1513, 154B/C, 1514D, 1515A/B, 1517A.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
PHILIPS IntelliVue Telemetry 3-ch	completo, one-piece		EE203SK-IV3	
PHILIPS IntelliVue Telemetry 5-ch	completo, one-piece		EE203SK-IV5	
PHILIPS StressVU.	completo, one-piece		EEHL20XL	
PHILIPS TC20 (no other TC models), C3i, Efcia ECG100.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
PHILIPS TC Set of 10 leadwires 4 mm Banana. Set 10 terminazioni banana 4 mm per Philips serie TC.	terminazioni, leadwires		EE103BK-TC	
PHILIPS TRIM & TOUCH Set of 10 leadwires 4 mm Banana. Set 10 terminazioni banana 4 mm.	terminazioni, leadwires		EE103BK-TRIM	
PHYSIO CONTROL Ready Link.	completo, one-piece	EE100MRI-200	EE100MRSI-200	EE100MRSI-350
PICKER VS-6.	completo, one-piece	EE100NRI-NV200		
POLYMED Alpha 6600, Alpha 8800.	completo, one-piece	EE100NRI-NV200		
PROGETTI conn. D-Sub 15 & MyECG, PG ECG.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
QMED Interp 1000.	completo, one-piece	EE100DRI-200	EE100DRSI-200	EE100DRSI-350
QUINTON BURDICK H12, H12+, X21, X12+	completo, one-piece		EEHL20XL	
ROZINN – SCOTT MEDICAL RZ151 2-channel / 5 DIN-Snap.	terminazioni, leadwires		EEHL05-60	
ROZINN – SCOTT MEDICAL RZ151 3-channel / 7 DIN-Snap, RZ152, RZ153+ 7D.	terminazioni, leadwires		EEHL07-60	
SAADAT Dena 350, 640, 650.	completo, one-piece	EE100ERI-200	EE100ERSI-200	

DESCRIZIONE / DESCRIPTION	TIPO / TYPE	TERMINAZIONI / TERMINATIONS		
		BANANA 4 MM	SNAP	SNAP XL (L = 4.26 m)
SCHILLER AT-1, AT-2, AT-2plus, AT-2 light, AT-4, AT-101, AT-104, AT-110. CS-6, CS-6/12, CS-12, AT- 5, AT- 6, AT-10, AT-10+, AT-60, AT-100. MS-2007, MS-2010, MS-2015. AT-10plus, AT-102, AT-104PC, CS-100, CS-200, CS200 excellence, CS200 touch etc.	completo, one-piece	EE100SRI-200	EE100SRSI-200	EE100SRSI-350
SCHILLER AT-3.	completo, one-piece	EE100NRI-NV200		
SCHWARZER Cardos 12, Light, CU12.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
SCHWARZER Cardos 6.	completo, one-piece	EE100CUBI-200	EE100CUSI-200	
SCHWARZER VS-6.	completo, one-piece	EE100NRI-NV200		
SECA Serie CT3000i, CT6i/6Pi, CT8000i/P/L/T. CT110, CT220. CT cardiopad; PT-160.	completo, one-piece	EE100SRI-200	EE100SRSI-200	EE100SRSI-350
SEIVA ECG Praktik, Pratik Portable, Cardiowriter, CardioTouch; Ergolog.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
SIEMENS Cardiostat 1, 11, 31, 31S.	completo, one-piece	EE100FRI-200	EE100FRSI-200	
SIEMENS Cardiostat 3T, 703, 710, 31, 31S; Cardirex 3T, 4, 6T, 62; Megacard Elema A/B, R, R/E; Cathcor, Cardina C, Miracard 1; Sirecard 440, 440S, 460; Mingodoc, Mingocard 3, 4, 7; Mingophon 3, 7, 310, 410, 420, 710, 720; Mingograph 4, 62, 81, 82, 740; Mingolog 7;	completo, one-piece	EE100CUBI-200	EE100CUSI-200	
SINOHERO SE101, SE505, SE508, SE600, SE1200, SE1200 Light.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
SONOSCAPE models conn. D-Sub 15.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
SONOSTAR SE-1, -1A, -3, -3B, -3C, -6, -6A, 12A, -12B, -12P, Ebox-12.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
SPACELABS Aria, Evo Ref. 25452-525.	terminazioni, leadwires		EEHL40	
SPACELABS CardioDirect CD12 USB, CD12-S USB; Cardiotrack 12, Cardiocollect.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
SPACELABS - BURDICK Holter 92513	cavo base, trunk cable		EEHL50	
SPACELABS – REYNOLDS Lifecard CF 3-ch: 7 DIN-Snap for trunk cable (trunk cable not included).	terminazioni, leadwires		EEHL07-60	
SPACELABS – DEL MAR Model 423, 423A,459, 461, 461A,463, 464, 483, 485.	terminazioni, leadwires		EEHL07-60	
SPACELABS – DEL MAR Model 419, 420,445, 447, 448,449, 453, 456.	terminazioni, leadwires		EEHL05-60	
SPACELABS - REYNOLDS Lifecard CF (3-ch x 3 leadwires)	completo, one-piece		EEHL30	
SPACELABS - REYNOLDS Lifecard CF (3-ch x 3 leadwires);	completo, one-piece		EEHL30/CLIP	

DESCRIZIONE / DESCRIPTION	TIPO / TYPE	TERMINAZIONI / TERMINATIONS		
		BANANA 4 MM	SNAP	SNAP XL (L = 4.26 m)
SPACELABS CardioExpress SL3, SL6, SL12.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
SPACELABS Voyager V6, V12.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
SPACELABS 90205, 90208, 92513.	terminazioni, leadwires		EEHL07-60	
SPACELABS Ultraview SL2400, 92513.	terminazioni, leadwires		EEHL05-60	
ST ELECTROMEDICINA 531, 531i, 563, 563i. ar600adv, ar600view, ar1200adv, ar1200view ...	completo, one-piece	EE100DRI-200	EE100DRSI-200	EE100DRSI-350
TEPA EKG Master.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
THOR Thoriox.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
TRISMED Cardicode 300 Type HT 312	terminazioni, leadwires		EEHL10-60	
TRISMED Cardipia all models, Stress-test models 406RS, 812RS.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
UNIVERSAL / UNIVERSALI Set of 5 pcs 60cm DIN-Snap ...	terminazioni, leadwires		EEHL05-60	
UNIVERSAL / UNIVERSALI Set of 7 pcs 60cm DIN-Snap ...	terminazioni, leadwires		EEHL07-60	
UNIVERSAL 5-LDS CABLE FOR VET & HUMAN / UNIVERSALE 5 TERM. PER VETERINARIA E UMANA.	completo, one-piece	EE100VET-200		
UTAS Ucard 100, Ucard 200.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
VALES & HILLS CV200. CV1200.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
WELCH ALLYN Cardio Perfect.	completo, one-piece	EE100CUBI-200	EE100CUSI-200	
WELCH ALLYN CP50, 150, 300.	completo, one-piece	EE100URI-200	EE100URSI-200	EE100URSI-350
WELCH ALLYN Serie AT. Serie SP. Serie CT, CP50, CP150, CP300.	completo, one-piece	EE100SRI-200	EE100SRSI-200	EE100SRSI-350
ZIMMER Cardio 1, Cardio 3	completo, one-piece	EE100CUBI-200	EE100CUSI-200	
ZIMMER Cardio 5, Cardiocompact.	completo, one-piece	EE100SRI-200	EE100SRSI-200	EE100SRSI-350
ZONCARE ZQ-1201, -1201G, -1203, -1203G, -1203S, -1206, -1212, IMAC730, T12.	completo, one-piece	EE100ERI-200	EE100ERSI-200	
ZYMED – PHILIPS 1100-009, 1100-010	terminazioni, leadwires		EEHL05-60	
ZYMED – PHILIPS 1100-012.	terminazioni, leadwires		EEHL07-60	

DESCRIZIONE / DESCRIPTION	TIPO / TYPE	TERMINAZIONI / TERMINATIONS		
		BANANA 4 MM	SNAP	SNAP XL (L = 4.26 m)