


## EE100GRI-200 PATIENT CABLE FOR RESTING ECG



NB the product is subject to changes without prior notice.

### Description

One-piece patient cable for Rest-ECG (12-leads diagnostic ECG). It is completely compatible with the original ones.

The plastic connector towards the device is conformal to the original. The terminations have 4mm Banana plugs.

Made for GE Marquette devices. The Compatibility List could be extended to further devices, according to the specific features.

Every trademark is propriety of the respective owner.

### Technical Features

<b>Length:</b>	✓ Main cable: 200 cm ✓ Leadwires: chest 70 cm, limbs 110 cm
<b>Terminations:</b>	10 pcs. 4mm Banana with isolated tips
<b>Yoke:</b>	compact, sealed
<b>Connector:</b>	D-sub 15 pin, made from plastic
<b>Colour code:</b>	IEC
<b>Resistors:</b>	10 k $\Omega$
<b>Protection against defibrillation:</b>	present
<b>Compatibility:</b>	GE Marquette MAC 400, 800, 1600, 3500, 500, 5500. GE Marquette Microsmart, MAC 500, 1000, 1100, 1200. GE style connector, not Hellige style.

