


EE100NRI-NV200 PATIENT CABLE FOR REST-ECG



NB the product is subject to changes without prior notice.

Description

One-piece patient cable for Rest-ECG (12-leads diagnostic ECG). It is completely compatible with the original ones. The plastic connector towards the device is conformal to the original. The terminations have 4mm Banana plugs..

The Compatibility List could be extended to further devices, according to the specific features.

Every trademark is propriety of the respective owner.

Technical Features

Length:	<ul style="list-style-type: none"> ✓ Main cable: 200 cm ✓ Leadwires: chest 70 cm, limbs 110 cm
Terminations:	10 pcs. 4mm Banana with isolated tips
Yoke:	compact, sealed
Connector:	D-sub 15 pin, made in plastic
Colour code:	IEC
Resistors:	4.7 kΩ
Defi protection:	present
Compatibility:	ASTRON Replay, Replay 3. BOSCH 501A, 501NV, 503, 503B, 503D, 503NV, 503V, 506, 506D, 601. BURDICK E350, E550, E560, EK10, Elite, Elite II . CAMBRIDGE VS-6, SC1000. CARDIETTE, ELETTRONICA TRENITINA every model of Micro, Microruler, Microtel. CARDIOLINE E1, E1R, E3, Delta300, Delta360, Delta1, Delta3. CARDIONICS other models with D-sub15 connector. COLSON DUPONT All models from the series Microruler. DR. LEE ECG-120A, 310A, 310B. FUKUDA DENSHI Autocardiner: FCP11, FX101. MCP302. Cardimax: FX2111, FCP2155, FX4010, FX7402. Alpha 500, Alpha 600, FCP 2201. FUKUDA ME Cardiosunny 1000, 6000A, 6600AX. HELDIGE EK51, Corina. HONEYWELL PHILIPS Cardioplan ER330, CR300, CR330. MEGOS Replay, Replay 3. MORTARA Eli100, 200. NIHON KOHDEN Cardiofax 6151, 6353, 8110, 8240, 8270, 8330, 8340, 8350, 8370, 8420, 8423, 8450, 8453, 8820, 8823, 8830. Cardiofax Q. PICKER VS-6. POLYMED Alpha 6600, Alpha 8800. SCHILLER AT-3. SCHWARZER VS-6 ...

