


EE100VET-200 CAVO PAZIENTE PER ECG A RIPOSO



NB il prodotto potrebbe differire da quanto illustrato. Il produttore si riserva il diritto di modificare il prodotto senza darne preavviso.

Descrizione

Cavo paziente intero per elettrocardiogramma a riposo a 5 fili, recanti 4 derivazioni periferiche e una precordiale. I cavi EE100VET-200 sono destinati all'uso in **Medicina Umana e Veterinaria**.

Il connettore lato apparecchio è realizzato in plastica ed è conforme all'originale. I terminali sono del tipo a banana 4 mm.

Realizzato come cavo universale per molteplici marche e modelli di elettrocardiografi.

La compatibilità può essere estesa ad altri apparecchi, in base alle caratteristiche specifiche dell'apparecchio in esame.

Tutti i marchi registrati o depositati sono proprietà dei rispettivi titolari.

Dati Tecnici

Lunghezza:	✓ Cavo base: 200 cm ✓ Terminazioni: precordiali 100 cm, periferiche 120 cm
Terminali:	5 terminali a banana 4 mm isolate R, L, F, N, V
Partitore:	pressofuso, compatto
Connettore:	D-sub 15 pin, in plastica
Codifica colori:	IEC
Resistenze:	10 kΩ
Protezione da defibrillazione:	presente
Compatibilità:	AMBISEA AV-9100, 9300, 9600, 9120. AMEDTEC ECGpro, CardioPart 12 Blue, CardioPart 12 USB. AVIONICS DEL MAR 910, 920, 930, 940, 942, 950, 960, 970, 980. ATES EBNEURO Easy ECG Pocket & Mobile. BIONET cardiocare 2000, cardiotouch 3000, cardio 7. BTL Serie BTL-08; CardioPoint-ECG C300, C600. CardioPoint-ERGO C300, C600. BOSCH 501A, 501NV, 503, 5038, 503D, 503NV, 503V, 506, 603D, 603ER, 606. BURDICK QUINTON Eclipse 800, Eclipse LE/LEi (> SN 1051), CareCenter MD. CAMEDIO 12B, 12H, 12S. CARDIOMAP 100H, 300H, 600H, 1200H. CARDIONICS CarTouch, Cardio Plug. CARDIOTEK serie EP-Tracer. CAREFUSION JAEGER Masterscreen ECG. COLIN YM406i. CORSCIENGE BT12 wireless. CUSTO.MED norm, cardio130. DAHIAN Araz 630. DR. SCHMIDT SmartScript. DR. LEE 310A, 310B. DR. VETTER pc-ekg. FIDELITY LP3000. FUKUDA ME Alpha 1000, 6600A, 6600AX. GEM-MED Serie Heart One +. GIMA Cardiopocket, Smart SE-1, -3, -6, -12 Express. Gimaplay 1. Cardiogima 1M, 3M, 6M, 12, 12M, Colour ECG. PC-ECG 33320, 33321; Acquisition

box series MGY. HELIGE EK53. HEWLETT PACKARD PageWriter 100, 200, 300 (M1770A, M1771A, M1772A). HORMANN VON BERG Bioset 9000. INNOMED Serie Heart Mirror & Heart Screen. Cardio PC/E. KENZ 108 & 110. LABTECH EC-12RT. M4MEDICAL M-trace, M-trace Mini, M-trace pc-ecg. MAQUAIRE CT-100, -200, -200i, -200CV. MEDEA MedeaLine 3, MedeaLine 6p, Kardio PCM-u, Excercise Test, Cardio Test. MEDIANA YM303i, 406i, 412i, 812i. MEDICAL ECONET Cardio-M, Cardio-M plus, Cardio-M PC, (instead of Quickels suction system / al posto del sistema a suzione). MEDIGATE ECG-1211-DG; MeCA303i, 406i, 412i, 812i. MEDISPO EM9530, 9535, 9540, 8033, 9550, 9560. MEDSET Serie Flashlight D-sub15 connector. MESA Cardiax pc-ecg. MS-WESTFALIA LunoCard S & HR and all models with D-sub15 connector. NEUROSOFT Serie Poly-Spectrum: -12/E, -8/EX, -8/E, -8,-Velo. PHILIPS HP PageWriter 10, 10i, Pagewriter 1513, 154B/C, 1514D, 1515A/B, 1517A. PHILIPS TC20 (no other TC models), C3i, Eficia ECG100. PROGETTI all models with D-sub15 connector & MyECG, PG ECG. SCHWARZER Cardos 12, Light, CU12. SEIVA ECG Praktik, Pratik Portable, Cardiowriter, CardioTouch; Ergolog. REYNOLDS SPACELABS CardioDirect CD12 USB, CD12-S USB; Cardiotrack 12, Cardicollect. SPACELABS Voyager V6, V12. THOR Thoriox. TEPA EKG Master. TRISMED Cardipia all models, stress-test models 406RS, 812RS. UTAS Ucard 100, Ucard 200. WELCH ALLYN CP50, 150, 300. Serie AT. Serie SP. Serie CT, Cardio Perfect . ASTRON 3/12b, 3/12i; Comp1x12, Comp3x12, Comp6x12, Comp3, Comp12, Comp60. CARDIETTE / ELETTRONICA TRENINA all versions of AR600,1200, 2100; Start 100P, Start 100H, Start 200, 200HV. DOC 101, 103; Autoruler 12/1, 12/3; Excel 103, 106, 106SCP; Excel 112, 212; Daedalus 346; Serie Daedalus View. CARDIOLINE all versions of AR600,1200, 2100; E1, E1R, E3; ECG100+, ECG200+. CARDIORAPID K11, K111, K112pocket, K130B, K130D, K131, K30P, K300P, K31, K312P, K33M, K36M, K360B, K360D, K3600LCD, K600P, K600M, K612. CARDIOSWIFT K112LCD, K36LCD & all models with D-sub 15 connector. COLSON DUPONT Cardi-3, E1, E3. Cardiopratique. All models from these series: AR, Cardipocket, Cardio, Excel, Autoruler, Microruler; Daedalus 346. ERGOLINE Autoruler: Autoscript 1/12, 3/12, 6/12 . Ergoscript EK301Z, EK1012, EK3003, EK3006, EK3012. FARUM all models with D-sub15 connector. MEGOS 3/12b, 3/12i; Comp1x12, Comp3x12, Comp6x12, Comp3, Comp12, Comp60. QMED Interp 1000. ST ELECTROMEDICINA 531, 531i, 563, 563i. ar600adv, ar600view, ar1200adv, ar1200view ...



CND Z12050380 Repertorio 212574/R.