


EE100VET-200 PATIENT CABLE FOR REST-ECG EE100VET-350 PATIENT CABLE FOR REST-ECG "XL"



NB the product is subject to changes without prior notice.

Description

EE100VET-200 and -350 are designed for Veterinary and Human Medicine. One-piece patient cable whose plastic connector towards the device is conformal to the original. The terminations have 4mm Banana plugs. The Compatibility List could be extended to further devices, according to the specific features.

Every trademark is propriety of the respective owner.

Tecnica Features

	EE100VET-200	EE100VET-350
Length:	<ul style="list-style-type: none"> ✓ Main cable: 200 cm ✓ Leadwires: chest 100 cm, limbs 120 cm 	<ul style="list-style-type: none"> ✓ Main cable: 350 cm ✓ Leadwires: chest 76 cm, limbs 126 cm
Terminations:	5 pcs. 4mm Banana with isolated tips	10 pcs. 4mm Banana, isolated tips
Yoke:	compact, moulded or sealed	
Connector:	D-sub 15 pin, made in plastic	
Colour code:	IEC	
Resistors:	10 kΩ	
Protection against defibrillation:	present	
Compatibility:	<ul style="list-style-type: none"> ✓ EE100VET-200: Universal protected cable with IEC leads 4mm Banana (R, L, F, N, Vx) and standard 15-pin connector. ✓ EE100VET-350: Universal protected cable with IEC leads 4mm Banana (12-Leads Standard) and standard 15-pin connector. 	

