

EEHL20xx patient cable for Holter, Stress-Testing & Monitoring



NB the product is subject to changes without prior notice.

Description

12-leads Holter or Stress-Testing patient cable. Unparalleled compatibility with the original cables. The connector is strong and makes a perfect match with the recorder. The main cable has progressive yoke with 12 Snap terminations. Made for Mortara, Cosmed, Philips and Datascope-Mindray devices. The Compatibility List could be extended to further devices, according to the specific features.

Every trademark is propriety of the respective owner.

Technical Features

Terminations:	10 Snap studs
Yoke:	progressive
Connector:	Mortara style.
Codifica colori:	IEC
Resistors:	not present
Protection against defibrillation:	not present
Packaging:	one piece each.
Compatibility:	COSMED Quark C12, T12. DATASCOPE Spectrum, Passport 2. MINDRAY Datascope Spectrum & Passport 2. MORTARA X-12, H-12, X-12+, H-12+ (no T12, T12+); Surveyor II / X-12 Transmitter. PHILIPS StressVU. QUINTON BURDICK H12, H12+, X21, X12+ ...

