

EEHL30 patient cable Lifecard CF EEHL30/CLIP patient cable for Lifecard CF, with belt clip



NB the product is subject to changes without prior notice

Description

3-leads Holter patient cable. Unparalleled compatibility with the original cables. The connector is strong and makes a perfect match with the recorder. The main cable has compact yoke with 3 Snap terminations.

Available with or without belt clip:

- EEHL30 corresponds to the original ref. 46-0480
- EEHL30 corresponds to the original ref. 46-0418

Every trademark is propriety of the respective owner.

Technical Features

Terminations:	3 Snap leads, L= 60 cm
Yoke:	molded
Connector:	- EEHL30: Lifecard. - EEHL30/CLIP: Lifecard with belt clip
Codifica colori:	IEC
Resistors:	not present
Protection against defibrillation:	not present
Packaging:	one piece each.
Compatibility:	ESAOTE Lifecard CF (3-ch x 3 leadwires); SPACELABS REYNOLDS Lifecard CF (3-ch x 3 leadwires)

