

EEHL50 patient cable Burdick Spacelabs 92513



NB the product is subject to changes without prior notice

Description

3-leads Holter patient cable for Burdick Spacelabs 92513.
Unparalleled compatibility with the original cables. The connector makes a perfect match with the recorder.
The trunk cable has a yoke for 7 DIN leadwires and it is available with the following color-coding: CH 1: Red-White, CH 2: Black-Brown, CH 3: Blue-Yellow/Orange, N: Green.
Can be used along with EEHL07-60.
The Compatibility List could be extended to further devices, according to the specific features.

Every trademark is propriety of the respective owner.

Technical Features

Leads:	3 bipolar leads plus Reference lead
Length:	60 cm
Yoke:	molded, for DIN leadwires
Connector:	proprietary
Color coding:	CH 1: Red-White, CH 2: Black-Brown, CH 3: Blue-Yellow/Orange, N: Green.
Resistors:	10k
Protection against defibrillation:	not present
Packaging:	one piece each.

