

EEHLxxx DIN-SNAP LEADWIRES



NB the product is subject to changes without prior notice.

Description

Set of leadwires for Holter, Telemetry and any kind of device with DIN sockets. Flexible and strong leadwire with Snap termination and DIN connector. Designed several Brand and models. The Compatibility List could be extended to further devices, according to the specific features.

Every trademark is propriety of the respective owner.

Technical Features

Terminations, connector:	Snap, DIN connector
Length:	60 cm
Packaging:	set of 3, 4, 5, 7 o 10 leadwires each, IEC Colours: <ul style="list-style-type: none"> ✓ EEHL03-60: Red, Yellow, Green ✓ EEHL04-60: White, Red, Orange, Brown ✓ EEHL05-60: Red, Black, White, Green, Yellow ✓ EEHL07-60: Red, Black, White, Green, Yellow, Blue, Brown ✓ EEHL10-60: Red, Black, Green, Yellow, Red, Yellow, Green, Brown, Black, Purple (Blue) ✓ Other combinations are available on request
Resistors:	not present
Defibrillation:	not protected against defibrillation discharge
Compatibility:	BMS Biomedical Software 5 DIN-to-snap leadwires, 7 DIN-to-snap leadwires. Braemar 5 DIN-to-snap leadwires, DL-400, DL-450, DL-700, DXP-1000. Cardiette - Elettronica Trentina DXP1000, Microscan 2000, Hr One, Hr One48. Cardioline Click Holter. Carewell CHS-12 Dynamic ECG. ContecTLC4000, TLC5000, TLC6000, TLC49804, TLC9803, PM85, PM60B. Datrix Intricon XR-300 and others with 5 or 7 DIN-to-snap leadwires. Del MarDigicorder 461, 483. DMS DM Systems, Meigaoyi DMS300-5, DMS 300-7, DMS 300-7P, DMS 300-7L, DMS 300-540, DMS300-8, DMS 300-8, DMS 300-12L. Esaote Evo, Aria, Lifecard CF 3-ch/7-leads (trunk cable not included). GE Marquette Hellige 3-channel recorders serie 8000, Seer, Seer MC; 2-channel recorders serie 7000, 8000, Seer, Seer XT, Seer MC. Gima ECG Holter 35130. Medilog (Oxford, Schiller) FD5. North-East Monitoring 5 or 7 DIN-to-snap leadwires. Oxford MedilogExcel, FD-4. Rozinn, Scott Medical RZ151 2-channel / 5 DIN-to-snap leadwires, RZ151 3-channel / 7 DIN-to-snap leadwires, RZ152, RZ153+ 7D. Spacelabs Ultraview SL2400, 92513, 90205, 90208, 92513, 90341 telemetry, Evo. Spacelabs, DelMar-Reynolds Aria, Lifecard CF 3-ch (trunk cable not included). Trismed Cardicode 300 Type HT 312. Zymed1100-009, 1100-010, 1100-012...

