



**"EUROFILM" Thermal Paper Roll for Video Printers**



**DESCRIPTION**

LUMED "EUROFILM" offers two models of high-resolution thermal paper roll for video printers (SONY UPP-110HD, SONY UPP-110HG and MITSUBISHI K65 HM-CE). They are used for printing high-resolution ultrasound scan images. Thanks to the optimal grey-scale response, our product solves the major problem of the other compatible films, because it can catch all the clinical information as the original does.

Each roll is individually wrapped in aluminum foil, each box contains 25 rolls.

**TECHINCAL FEATURES**

REFERENCE	BMSY013	BMSY014
DESCRIPTION	Thermal Paper Roll for Video Printers	
COMPATIBILITY	UPP-110HD and K65 HM-CE	UPP-110HG and K65 HM-CE
DIMENSIONS	110 x 20mt	110 x 18mt
SALE PACKAGING	25 rolls in a box	25 rolls in a box
PACKAGING MATERIAL	Aluminum Foil	
COMPLIANCE TO STANDARDS	<p><b>Directive 93/42/EEC and amendments. Class I medical device</b>  <b>Italian Medical Device Data standard Registration:</b>  <b>56383/R; CND: Z 1302; GMDN: 15639; CE: 0476</b>                      *Non-Active Medical Device*</p>	
MANUFACTURER	<p><b>LUMED srl</b>                      Via Staffora, 18/9 - 20090 Opera (MI) - ITALY</p>	
STORAGE	<p>Temperature: &lt; 40°C – Relative Humidity: &lt; 65%                      Do not expose the film to reactive gases, sun light, UV sources etc.</p>	